

Protocol of BAC Fingerprinting with SNaPshot™ Kit

1. **Re-suspension:** BAC DNAs isolated by Qiagen R.E.A.L Prep 96 Plasmid kit were re-suspended in 42 μ l of ddH₂O in deep well block, vortex gently, and keep at room temperature overnight or at 4 °C over weekend..
2. **Restriction enzyme digestion:** Add 8.0 μ l of the restriction enzyme cocktail (see below); incubate at 37 °C, for 3 hrs or longer (overnight is OK).

Enzyme Cocktail (1x) (NEB enzymes):

<i>Bam</i> HI	2.0 units (0.10 μ l)
<i>Eco</i> RI	2.0 units (0.10 μ l)
<i>Xba</i> I	2.0 units (0.10 μ l)
<i>Xho</i> I	2.0 units (0.10 μ l)
<i>Hae</i> III	2.0 units (0.20 μ l)
NEBuffer 2	5.0 μ l
100X BSA	0.5 μ l
RNase A (0.5 μ g/ μ l, DNase free)	1.0 μ l
β -Mercaptoethanol (1%)	1.0 μ l

3. Transfer 50.0 μ l of the digested DNAs into 96-PCR plate.
4. **Labeling:** Add 10.0 μ l of SNaPshot labeling cocktail (see below), briefly spin down; incubate at 65 °C for 60'.

Labeling Cocktail (1x):

SNaPshot Multiplex Ready Reaction Mix (from ABI)	0.3 μ l
NEBuffer 2	2.0 μ l
100 mM Tris (pH = 9.0)	2.5 μ l
ddH ₂ O	5.2 μ l

5. **Precipitation:** Add 5.0 μ l of 2.5M Sodium Acetate, 100 μ l of pre-chilled ethanol (95%), and place at -80°C for 10-15'. Spin at 4200 rpm for 30'; Wash with 70% ethanol and spin at 3500 rpm for 10'; spin upside down on paper towel at 300 rpm for 2'; air dry for 5' or longer.
6. **Re-suspension:** Re-suspend the pellet with mixture of 9.85 μ l of Hi-Di formamide and 0.15 μ l of GS1200Liz Size Standard, then vortex gently.
7. **ABI 3730XL:** Denature at 95 °C for 5', and place on ice until ready to load on the ABI 3730.